


PRODUCT BROCHURE

KERN 346-000-000

Date: 17.04.2025

ENGLISH



Individual weight container for rectangular weights

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- Individual weight containers for testing high-load floor scales, pallet scales, pallet truck scales, crane scales, etc. This can also be used for storing the weights.
- This means the weight container and the weights can be placed on the balance in one go, saving you time and money
- The weight container can be calibrated to OIML accuracy classes M1 – M3
- On request, KERN will make you a “tailor-made” weight carrier to your specifications

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	346-000-000
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1224995
Customs tariff number	84239090

Dimensions packaging (W×D×H)	1000×500×500 mm
------------------------------	-----------------

Shipping method	Parcel service
-----------------	----------------

Delivery time	1 d
---------------	-----

Optional Services

DAkkS Calibration	962-600
-------------------	---------

FUNCTIONS

Standard



Technical data - Construction

Dimensions	1000×500×500 mm
------------	-----------------

Material	stainless steel
----------	-----------------

Technical data - Measuring system

OIML class	M1
------------	----

OIML conform	yes
--------------	-----

Technical data - Packing & Shipping

SERVICES

Model	Description
962-646	DAkkS Calibration KERN 962-646
962-600	DAkkS Calibration KERN 962-600
962-600R	DAkkS Recalibration KERN 962-600R

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium
346-06	5 kg	250 mg	346-060-600
346-07	10 kg	500 mg	346-070-600
346-08	20 kg	1 g	346-080-600
346-09	50 kg	2,5 g	346-090-600
346-022-005	20 kg	1 g	-
346-042-008	40 kg	1,5 g	-
346-052-010	50 kg	2,5 g	-
346-055-004	50 kg	2,5 g	-
346-055-009	50 kg	2,5 g	-
346-065-009	60 kg	3 g	-
346-86	5 kg	250 mg	346-060-600
346-87	10 kg	500 mg	346-070-600
346-77	10 kg	500 mg	346-070-600
346-76	5 kg	250 mg	346-060-600
346-79-2023a	50 kg	2,5 g	-
346-88	20 kg	1 g	346-080-600
346-78	20 kg	1 g	346-080-600
346-89	50 kg	2,5 g	-
346-81	100 kg	5 g	-
346-82	200 kg	10 g	-
346-83	500 kg	25 g	-
346-84	1000 kg	50 g	-
346-79	50 kg	2,5 g	346-090-600