

KERN MWA 300K-1PM

Wheelchair platform scale with two integrated access ramps for easy access – with approval for professional medical use in medical diagnostics, verification optional



Article codes, product groups & Status information

Model article-code	MWA 300K-1PM
Model-Serie	MWA
Assortment	Scales
Product Type	Medical products
Product group	Wheelchair platform scales
Product usage type	Main article
EAN code (Model)	4045761193048
Predecessor Model	MWA 300K-1PM-2.017a
Planned start date	2017-2HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary

Weighing capacity	300 kg
Readability	0,1 kg
Verification value (e)	0,1 kg

Technical data - Construction

Material weighing plate	Steel, lacquered
Bubble Level	Yes
Weighing surface (d)	880×840 mm
Dimensions of weighing plate (W×D)	880×840 mm
Dimension completely mounted (W×D×H)	1.180×1.130×1.180 mm

Material casing	Steel, lacquered
IP protection load cell	IP65

Technical data - Display

Display height	25 mm
----------------	-------

Technical data - Power supply

Battery	6×15 V AA
Operating time	20 h
Mains adapter/adaptor included	CH;EURO;UK
Charging time	12 h
Input Voltage	100 V - 240 V, 50 / 60 Hz
Operating time (Backlight off)	40 h
Rechargeable battery	Rchrg. battery optional - intern
Interfaces	RS-232 standard

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 0,1 kg
Readability	0,1 kg
Recommended adjusting weight	300 kg (M1)
Repeatability	0,1 kg
Weighing units	Kg
Weighing system	Strain gauge
Max. divergency at eccentric load (1/3 max)	0,1 kg

Technical data - Environmental conditions

KERN MWA 300K-1PM

Wheelchair platform scale with two integrated access ramps for easy access – with approval for professional medical use in medical diagnostics, verification optional



Maximum humidity	80 %
------------------	------

Maximum operating temperature	40 °C
-------------------------------	-------

Minimum ambient temperature	10 °C
-----------------------------	-------

Minimum humidity	5 %
------------------	-----

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	1.255×1.100×285 mm
------------------------------	--------------------

Gross weight	90 kg
--------------	-------

Net weight	80 kg
------------	-------

Shipping method	Freight
-----------------	---------

Packaging component - by weight - cardboard	7 kg
---	------

Packaging component - by weight - styrofoam	7 kg
---	------

Delivery time	2 d
---------------	-----

Technical data - Verification

Verification Value E	0,1 kg
----------------------	--------

Minimum load (Min)	2 kg
--------------------	------

Verification by KERN possible	Yes
-------------------------------	-----

Verification class	III
--------------------	-----

MPG class	1m
-----------	----

CND code	V0399
----------	-------

GMDN code	15053
-----------	-------

Optional Services

DAkkS Calibration	963-129
-------------------	---------

Verification	965-129
--------------	---------

Functions

Standard



Option



Factory

