

## PRODUCT BROCHURE

---

# KERN 347-116

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M1, 1 kg, Hook weight, stainless steel fine turned

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

- Individual weight, OIML Class M1, 1 kg, Hook weight, stainless steelfine turned (OIML)

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 347-116                        |
| Model-Serie           | 347-116                        |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML M1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761256248                  |
| Predecessor Model     | 347-516                        |
| Customs tariff number | 84239090                       |
| Country of Origin     | DE                             |

### Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 1 kg        |
| Tolerance (OIML) | 50 mg       |
| Case plastic     | 347-120-400 |

### Technical data - Construction

|          |                             |
|----------|-----------------------------|
| Material | stainless steel fine turned |
|----------|-----------------------------|

### Technical data - Measuring system

|                              |          |
|------------------------------|----------|
| Nominal value                | 1 kg     |
| OIML class                   | M1       |
| OIML conform                 | yes      |
| Tolerance (OIML)             | 50 mg    |
| Density (kg/m <sup>3</sup> ) | 8000     |
| Volume (cm <sup>3</sup> )    | 119,0476 |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 48×48×123 mm   |
| Gross weight                               | 1.00 kg        |
| Net weight                                 | 1.00 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

### Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-641 |
|-------------------|---------|

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



|        |
|--------|
| Option |
|--------|



## ACCESSORIES

| Model                              | Description                         |
|------------------------------------|-------------------------------------|
| <a href="#"><u>317-280</u></a>     | Gloves KERN 317-280                 |
| <a href="#"><u>347-120-400</u></a> | Plastic Weight Box KERN 347-120-400 |

## SERVICES

| Model                           | Description                       |
|---------------------------------|-----------------------------------|
| <a href="#"><u>962-641</u></a>  | DAkkS Calibration KERN 962-641    |
| <a href="#"><u>962-641R</u></a> | DAkkS Recalibration KERN 962-641R |