

KERN 336-38

Individual weight, OIML Class F2, 20 kg, Block, stainless steel



Article codes, product groups & Status information

Model article-code	336-38
Model-Serie	336
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761232822
Planned start date	2020-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	20 kg
Tolerance (OIML)	300 mg
Case aluminium	346-080-600
Case plastic	Not available
Case wood	Not available

Technical data - Construction

Dimensions	233×115×147 mm
Material	Stainless steel

Technical data - Measuring system

Nominal value	20 kg
OIML class	F2
OIML conform	Yes
Tolerance (OIML)	300 mg

Functions

Standard

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	233×115×145 mm
Gross weight	20 kg
Net weight	20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-445
-------------------	---------

KERN 336-38

Individual weight, OIML Class F2, 20 kg, Block, stainless steel

