

KERN 347-115-100

Supporting bar, OIML Class M1, 1 kg, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-115-100
Model-Serie	347-115-100
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761256378
Predecessor Model	347-515-100
Planned start date	2022-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 kg
Tolerance (OIML)	50 mg
Case plastic	347-515-400

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	1 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	50 mg
Density (kg/m ³)	8000
Volume (cm ³)	119,0476

Functions

Standard

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	60×60×991 mm
Gross weight	1.00 kg
Net weight	1.00 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-641
-------------------	---------

KERN 347-115-100

Supporting bar, OIML Class M1, 1 kg, stainless steel fine turned

