

KERN 366-76

Individual weight, OIML Class M3, 5 kg, Block Eco, cast iron lacquered



Article codes, product groups & Status information

Model article-code	366-76
Model-Serie	366
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761241824
Successor Model	366-76-2025e
Planned start date	2020-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	2,5 g
Case aluminium	346-060-600

Technical data - Construction

Dimensions	150×75×84 mm
Material	Cast iron lacquered

Technical data - Measuring system

Nominal value	5 kg
OIML class	M3
OIML conform	Yes
Tolerance (OIML)	2,5 g
Density (kg/m ³)	7100

Functions

Standard

Volume (cm³) 704,2254

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×75×84 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-643
-------------------	---------

KERN 366-76

Individual weight, OIML Class M3, 5 kg, Block Eco, cast iron lacquered

