

KERN 346-86-2024a

Individual weight, OIML Class M1, 5 kg, Block, cast iron lacquered



Article codes, product groups & Status information

Model article-code	346-86-2024a
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761241626
Predecessor Model	346-86-2020a
Successor Model	346-86; 346-86-2024e2
Planned start date	2020-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	250 mg
Case aluminium	346-060-600

Technical data - Construction

Dimensions (W×D×H)	150×80×90 mm
Material	Cast iron lacquered

Technical data - Measuring system

Nominal value	5 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	250 mg

Functions

Standard

Density (kg/m ³)	7100
Volume (cm ³)	704,2254

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	170×95×100 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-643
-------------------	---------

KERN 346-86-2024a

Individual weight, OIML Class M1, 5 kg, Block, cast iron lacquered

