

KERN 308-36

Milligram weights, E1, 50 mg, Wire, stainless steel polished



Article codes, product groups & Status information

Model article-code	308-36
Model-Serie	308
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761481015
Predecessor Model	308-36-2.024a
Planned start date	2025-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	50 mg
Tolerance (OIML)	0,004 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	12×12×0,5 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	50 mg
OIML class	E1

Functions

Standard

OIML conform	Yes
Tolerance (OIML)	0,004 mg
Density (kg/m ³)	8.000
Volume (cm ³)	0,0063

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	12×12×3 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-256
-------------------	---------

KERN 308-36

Milligram weights, E1, 50 mg, Wire, stainless steel polished



Option

