

## PRODUCT BROCHURE

---

# OPTICS OBB-A1110

Date: 16.03.2025

ENGLISH



Achromatic objective lens, 20 x /0,4 (sprung) W.D. (1,75 mm)

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)

## DESCRIPTION

Achromatic objective lens, 20 x /0,4 (sprung) W.D. (1,75 mm) OBB-A1110

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                 |
|-----------------------|---------------------------------|
| Model article-code    | OBB-A1110                       |
| Model-Serie           | OBB-AO                          |
| Assortment            | Optics                          |
| Product Type          | Microscopes                     |
| Product group         | Transmitted light microscopes   |
| Product usage type    | Spare part;Accessories/Services |
| EAN code (Model)      | 4045761164499                   |
| Customs tariff number | 90119000                        |

### FUNCTIONS

Standard



ENGLISH

### Technical data - Environmental conditions

|                  |       |
|------------------|-------|
| Storage Temp Min | -5 °C |
| Storage Temp Max | 40 °C |

### Technical data - Packing & Shipping

|                                             |                |
|---------------------------------------------|----------------|
| Dimensions packaging (W×D×H)                | 65×40×40 mm    |
| Gross weight                                | 0.60 kg        |
| Net weight                                  | 0.30 kg        |
| Shipping method                             | Parcel service |
| Packaging component - by weight - cardboard | 300 g          |
| Delivery time                               | 1 d            |