

# KERN HBD 100-0

Compact handheld durometer with drag indicator

**KERN**<sup>®</sup>



## Article codes, product groups & Status information

Model article-code	HBD 100-0
Model-Serie	HB
Assortment	Sauter
Product Type	Hardness measurement
Product group	Hardness testers (Shore)
Product usage type	Main article
EAN code (Model)	4045761140189
Predecessor Model	HBD 100-0-2.019a
Planned start date	2012-1HJ
Customs tariff number	90248000
Country of Origin	CN

## Technical data - Primary

Hardness scales	Shore D
-----------------	---------

## Technical data - Construction

Dimension housing (W×D×H)	115×60×25 mm
Material housing	Metal

## Technical data - Measuring system

Test force Hardness measurement	10 N
Readability (Shore)	1
Tolerance (% from Max.)	3%
Unit	HD

## Technical data - Environmental conditions

## Functions

Maximum operating temperature	40 °C
Minimum operating temperature	0 °C
Maximum storage temperature	40 °C
Minimum storage temperature	-10 °C

## Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×95×50 mm
Gross weight	0,50 kg
Net weight	0,20 kg
Shipping method	Parcel service
UN number	No dangerous goods
Delivery time	1 d

# KERN HBD 100-0

Compact handheld durometer with drag indicator



Standard

