

KERN MSF 200

When every centimetre counts – mechanical height rod approved as a medical device for professional use in medical diagnostics



Article codes, product groups & Status information

Model article-code	MSF 200
Model-Serie	MSF
Assortment	KERN
Product Type	Medical products
Product group	Body size measurement
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761192843
Predecessor Model	MSF 200-2016a
Planned start date	2017-1HJ
Customs tariff number	84239010
CE Label	Yes

Technical data - Primary

Measurement range	210 cm
Readout	0.1 cm

Technical data - Construction

Dimensions (W×D×H)	60×60×930 mm
--------------------	--------------

Technical data - Measuring system

Measurement range	205 cm
Readout	0.1 cm
Min Length	60 cm
Device type	Mechanical
Units	Cm;inch/ft
Tolerance Length	10 mm

Functions

Standard

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	1015×70×45 mm
Gross weight	1.2 kg
Net weight	0.80 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	108.5 g
Packaging component - by weight - plastics	108.5 g
Delivery time	1 d

Technical data - Verification

Medical approval 93/42/EWG	Yes
MPG class	1m
GMDN code	37001

KERN MSF 200

When every centimetre counts – mechanical height rod approved as a medical device
for professional use in medical diagnostics

