

LED gooseneck illumination OZB-A4516



Function

1. On/Off-Button
2. Dimmer light intensity



Advantages at a glance

With **the OZB-A4516** 20 W LED gooseneck illuminator with focusable light cone, you can adjust your illumination individually. A punctual or a scattered beam enables you to optimally illuminate the specimen.



Safety informations

The light from the gooseneck illumination is very bright. Therefore, always switch on the illumination before looking through the eyepieces. Avoid looking directly into the LEDs.



Technical details

Illuminance: 20W-LED
 Colour temperature: 6400 K
 Length: 540 mm
 Input voltage : AC 100-240V 50/60HZ



Techn. service +49 7433 9933 777