

## PRODUCT BROCHURE

---

# OPTICS OBB-A1655

Date: 11.03.2026

ENGLISH



Image not  
available

Infinity planachromatic objective 2x /0,06 W.D. (7,2 mm)

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	OBB-A1655
Model-Serie	N/A
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761495166
Customs tariff number	90119000
Country of Origin	CN

Standard



ENGLISH

### Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

### Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	72×42×42 mm
Gross weight	0.15 kg
Net weight	0.15 kg
Shipping method	Parcel service
Delivery time	1 d

## FUNCTIONS