

KERN 308-33

Milligram weights, E1, 5 mg, Wire, stainless steel polished



Article codes, product groups & Status information

Model article-code	308-33
Model-Serie	308
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761480995
Predecessor Model	308-33-2.024a
Planned start date	2025-2HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	5 mg
Tolerance (OIML)	0,003 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	11×11×0,3 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	5 mg
OIML class	E1
OIML conform	Yes

Tolerance (OIML) 0,003 mg

Density (kg/m³) 8.000

Volume (cm³) 0,0006

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	11×11×3 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-253
-------------------	---------