

## PRODUCT BROCHURE

---

# KERN PCB-A01

Date: 16.03.2025

ENGLISH



Internal rechargeable battery pack, operating time up to 48 h, charging time approx. 8 h, AUTO-OFF function to preserve the rechargeable battery pack

---

## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	PCB-A01
Model-Serie	PCB
Assortment	KERN
Product Type	Laboratory balances
Product group	Accessories (printer, etc.)
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761085657
Predecessor Model	PCB-A06
Customs tariff number	85075000
CE Label	yes

### Technical data - Power supply

Mains adapter/adapter included	CH;EURO;UK;US
Mains adapter	Battery & Power supply

### Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	70×100×100 mm
Gross weight	0.40 kg
Net weight	0.25 kg
Shipping method	Parcel service
UN number	3496
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	82.92 g

Packaging component - by weight - plastics 82.92 g

Packaging component - by weight - styrofoam 0 g

Packaging component - by weight - glass 0 g

Packaging component - by weight - wood 0 g

Packaging component - by weight - foam 0 g

Packaging component - by weight - metal 0 g

Delivery time 1 d

## FUNCTIONS

Standard

