

KERN 347-137

Individual weight, OIML Class M1, 5 kg, Slot weight (Type B), stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-137
Model-Serie	347-137
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366831
Planned start date	2024-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	250 mg

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	5 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	250 mg
Density (kg/m ³)	8.000
Volume (cm ³)	595,2381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	178×178×34 mm
Gross weight	5,0 kg
Net weight	5,0 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-643
-------------------	---------

Functions

KERN 347-137

Individual weight, OIML Class M1, 5 kg, Slot weight (Type B), stainless steel fine turned



Standard



Option

