

KERN FOB 500-1S

Compact stainless steel balance for weighing, checking and portioning in the smallest of spaces



Article codes, product groups & Status information

Model article-code	FOB 500-1S
Model-Serie	FOB
Assortment	KERN
Product Type	Food scales
Product group	Bench scales
Product usage type	Main article
EAN code (Model)	4045761135734
Planned start date	2012-2HJ
Customs tariff number	84238129
CE Label	Yes

Technical data - Primary

Weighing capacity	500 g
Readability	0.1 g

Technical data - Construction

Material weighing plate	Stainless steel
Weighing surface (d)	125×155 mm
Dimensions of weighing plate (W×D)	125×155 mm
Material casing	Stainless steel

Technical data - Display

Display height	20 mm
----------------	-------

Technical data - Power supply

Battery	9 V Block
Operating time	20 h

Input Voltage	230 V AC, 50 Hz
---------------	-----------------

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 0.2 g
Readability	0.1 g
Recommended adjusting weight	500 g (M1)
Repeatability	0.1 g
Weighing units	Dwt;g;lb;oz;ozt
Number of weighing ranges	1
Weighing system	Strain gauge

Technical data - Environmental conditions

Maximum humidity	95 %
Maximum operating temperature	40 °C
Minimum ambient temperature	10 °C
Minimum humidity	25 %

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	220×190×75 mm
Gross weight	0.85 kg
Net weight	0.70 kg
Shipping method	Parcel service
UN number	3496
Packaging component - by weight - paper	100 g

KERN FOB 500-1S

Compact stainless steel balance for weighing, checking and portioning in the smallest of spaces



Packaging component - by weight - cardboard	82 g
---	------

Packaging component - by weight - plastics	82 g
--	------

Packaging component - by weight - styrofoam	0 g
---	-----

Packaging component - by weight - glass	0 g
---	-----

Packaging component - by weight - wood	0 g
--	-----

Packaging component - by weight - foam	0 g
--	-----

Packaging component - by weight - metal	0 g
---	-----

Delivery time	1 d
---------------	-----

Optional Services

DAkkS Calibration	963-127
-------------------	---------

Functions

Standard



Option

