

KERN 356-77

Individual weight, OIML Class M2, 10 kg, Block Eco, cast iron lacquered



Article codes, product groups & Status information

Model article-code	356-77
Model-Serie	356
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761241794
Planned start date	2020-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	1,6 g
Case aluminium	346-070-600

Technical data - Construction

Dimensions	190×95×109 mm
Material	Cast iron lacquered

Technical data - Measuring system

Nominal value	10 kg
OIML class	M2
OIML conform	Yes
Tolerance (OIML)	1,6 g
Density (kg/m ³)	7100
Volume (cm ³)	1408,4507

Functions

Standard

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	190×95×109 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-644
-------------------	---------

KERN 356-77

Individual weight, OIML Class M2, 10 kg, Block Eco, cast iron lacquered

