

## PRODUCT BROCHURE

---

# KERN 323-036

Date: 10.03.2026

ENGLISH



Weight set, F1, 1 mg - 100 g, Knob, stainless steel polished, in aluminium case

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                      |
|-----------------------|----------------------|
| Model article-code    | 323-036              |
| Model-Serie           | 323                  |
| Assortment            | KERN                 |
| Product Type          | Test weights         |
| Product group         | OIML F1 Test weights |
| Product usage type    | Main article         |
| EAN code (Model)      | 404576 12238 13      |
| Customs tariff number | 84239090             |
| Country of Origin     | DE                   |

### Technical data - Primary

|               |              |
|---------------|--------------|
| Nominal value | 1 mg - 100 g |
|---------------|--------------|

Case aluminium 313-030-600

Weight set case in aluminium case

### Technical data - Construction

|               |                          |
|---------------|--------------------------|
| Material      | stainless steel polished |
| Is weight set | yes                      |
| Type of case  | in aluminium case        |

### Technical data - Measuring system

|               |              |
|---------------|--------------|
| Nominal value | 1 mg - 100 g |
| OIML class    | F1           |
| OIML conform  | yes          |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 165×125×75 mm  |
| Gross weight                               | 0.90 kg        |
| Net weight                                 | 0.75 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 170 g          |
| Delivery time                              | 1 d            |

## Technical data - Verification

|                       |     |
|-----------------------|-----|
| Verification possible | yes |
|-----------------------|-----|

## Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-402 |
|-------------------|---------|

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



|        |
|--------|
| Option |
|--------|



## SERVICES

| Model                  | Description                                          |
|------------------------|------------------------------------------------------|
| <b><u>962-402</u></b>  | DAkkS Calibration KERN 962-402                       |
| <b><u>962-402R</u></b> | DAkkS Recalibration KERN 962-402R                    |
| <b><u>921-402</u></b>  | Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-402 |
| <b><u>920-703</u></b>  | Marking (Last digits of serial number) KERN 920-703  |

**SPARE PARTS**

| Model                              | Description                               |
|------------------------------------|-------------------------------------------|
| <u><a href="#">313-030-600</a></u> | Aluminium Protective Box KERN 313-030-600 |