

KERN ODC 874

Specialists in microscopy for measurement, counting, documentation, archiving and image processing



Article codes, product groups & Status information

Model article-code	ODC 874
Model-Serie	ODC-87
Assortment	Optics
Product Type	Microscopes
Product group	Microskope cameras
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761190528
Planned start date	2016-1HJ
Customs tariff number	85258900
CE Label	Yes

Technical data - Primary

Camera type	USB2 Eyepiece camera - 23.2 mm
Camera resolution	3,1 MP
Sensor type	Aptina CMOS
Sensor size	1/2,5"

Technical data - Construction

Color capability	Colour
Camera mounting	Plug-in connection
Interfaces	USB 2.0
Operating system	Win 10;Win 11;Win 8
Applicable software	Microscope VIS 2.0 Pro (OXM 902)
Dimensions (W×D×H)	40×27×27 mm

Technical data - Power supply

Protection Class	Class III - Low voltage
Cable attachment	removable

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	145×95×60 mm
Gross weight	0.50 kg
Net weight	0.45 kg
Shipping method	Parcel service
Packaging component - by weight - paper	72.1 mg
Packaging component - by weight - cardboard	34.5359 g
Packaging component - by weight - plastics	504.7 mg
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - foam	36.9873 g
Delivery time	1 d

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Verification

WEEE	Yes
------	-----

Functions

KERN ODC 874

Specialists in microscopy for measurement, counting, documentation, archiving and image processing



Standard

