


PRODUCT BROCHURE

KERN KUM-04

Date: 16.03.2025

ENGLISH



Internal data interface Ethernet

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- The KERN Universal Modules (KUM) enable the exchange of data and control commands between balance and network or PC, mobile end unit, printer etc. using several interfaces: RS-232, USB device, Ethernet, Bluetooth and WiFi. In this way we ensure the perfect integration of your KERN balance into production or processing chains, which gives the benefit of a simple, consistent workflow
- The KUM interface modules for RS-232, USB, Ethernet, Bluetooth and WiFi are plugged straight onto the board to suit individual choice and this means that the IP protection of the housing remains guaranteed

SPECIFICATIONS

Article codes, product groups & Status information

| | |
|-----------------------|------------------------------|
| Model article-code | KUM-04 |
| Model-Serie | KUM |
| Assortment | KERN |
| Product Type | General weighing accessories |
| Product group | Accessories (printer, etc.) |
| Product usage type | Accessories/DL |
| EAN code (Model) | 404576 1476752 |
| Customs tariff number | 90160090 |

Delivery time 1 d

FUNCTIONS

Standard



Technical data - Packing & Shipping

| | |
|---|----------------|
| Dimensions packaging (W×D×H) | 225×145×40 mm |
| Gross weight | 0.50 kg |
| Net weight | 0.30 kg |
| Shipping method | Parcel service |
| Packaging component - by weight - cardboard | 100 g |
| Packaging component - by weight - styrofoam | 100 g |