

PRODUCT BROCHURE

OPTICS OZB-A45 16

Date: 30.04.2025

ENGLISH



Double gooseneck illumination for best illumination of your sample

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



ENGLISH

DESCRIPTION

- With the 20 W-LED gooseneck illumination unit with focusable light beam, you can adjust the illumination to suit your needs. Spot or scattered radiation means that you can achieve the very best illumination of your sample

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OZB-A45 16
Model-Serie	OZB-IF
Assortment	Optics
Product Type	Microscopes
Product group	External Illumination

Product usage type	Main article
EAN code (Model)	404576 1217720
Predecessor Model	OZB-A45 12
Customs tariff number	94054239

Technical data - Primary

Illumination type	LED
Illumination intensity	20W

Technical data - Construction

Color tempertaure	6400 K
Dimmable	yes
Fibre length	540 mm
Dimensions (W×D×H)	100×200×150 mm

Technical data - Power supply

Protection Class	Class III - Low voltage
------------------	-------------------------

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	300×250×250 mm
Gross weight	4 kg
Net weight	3.8 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	207.576 g
Packaging component - by weight - plastics	1.41932 g
Packaging component - by weight - styrofoam	145.835 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Verification

WEEE	yes
------	-----

FUNCTIONS

Standard



SPARE PARTS

Model	Description
RKFN-5001	Power socket RKFN-5001
OZB-A4535	Replacement focus lenses (pair) OZB-A4535
OZB-A4581	Power supply KERN OZB-A4581
OZB-A4582	Power supply KERN OZB-A4582
OZB-A4534	Replacement light guide OZB-A4534
OZB-A4536	Replacement LED 20 W KERN OZB-A4536
RCFS-5001	Switch KERN RCFS-5001