

KERN 326-054

Weight set, F1, 1 g - 500 g, Eco Shape, stainless steel polished, in plastic case



Article codes, product groups & Status information

Model article-code	326-054
Model-Serie	326
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	4045761223950
Predecessor Model	322-054
Planned start date	2019-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 g - 500 g
Case plastic	313-052-400
Weight set case	In plastic case

Technical data - Construction

Material	Stainless steel polished
Is weight set	Yes
Type of case	In plastic case

Technical data - Measuring system

Nominal value	1 g - 500 g
OIML class	F1
OIML conform	Yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	250×210×75 mm
Gross weight	1.6 kg
Net weight	0.40 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	1080 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-418
-------------------	---------

Functions

Standard

KERN 326-054

Weight set, F1, 1 g - 500 g, Eco Shape, stainless steel polished, in plastic case

