

KERN 366-93

Individual weight, OIML Class M3, 500 g, Knob, cast iron lacquered



Article codes, product groups & Status information

Model article-code	366-93
Model-Serie	366
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761012080
Successor Model	366-93-2025e
Planned start date	2004-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	500 g
Tolerance (OIML)	250 mg

Technical data - Construction

Material	Cast iron lacquered
----------	---------------------

Technical data - Measuring system

Nominal value	500 g
OIML class	M3
OIML conform	Yes
Tolerance (OIML)	250 mg
Density (kg/m ³)	7100
Volume (cm ³)	70,4225

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	46×46×68 mm
Gross weight	0.50 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	No
-----------------------	----

Optional Services

DAkKS Calibration	962-639
-------------------	---------

Functions

KERN 366-93

Individual weight, OIML Class M3, 500 g, Knob, cast iron lacquered



Standard

