

KERN 326-066

Weight set, F1, 1 g - 1 kg, Eco Shape, stainless steel polished, in aluminium case



Article codes, product groups & Status information

Model article-code	326-066
Model-Serie	326
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	4045761224100
Planned start date	2019-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 g - 1 kg
Case aluminium	314-060-600
Weight set case	In aluminium case

Technical data - Construction

Material	Stainless steel polished
Is weight set	Yes
Type of case	In aluminium case

Technical data - Measuring system

Nominal value	1 g - 1 kg
OIML class	F1
OIML conform	Yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	230×140×110 mm
------------------------------	----------------

Functions

Gross weight	3.6 kg
Net weight	0.85 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	2690 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-419
-------------------	---------

KERN 326-066

Weight set, F1, 1 g - 1 kg, Eco Shape, stainless steel polished, in aluminium case



Standard

