

PRODUCT BROCHURE

KERN 354-046

Date: 10.03.2026

ENGLISH



Weight set, M2, 1 g - 200 g, Knob, stainless steel fine turned, in aluminium case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	354-046
Model-Serie	354
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main article
EAN code (Model)	404576 12923 14
Predecessor Model	354-446
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 200 g
Case aluminium	314-040-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 g - 200 g
OIML class	M2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	190×130×80 mm
Gross weight	1.2 kg
Net weight	0.90 kg
Shipping method	Parcel service

Packaging component - 200 g
by weight - plastics

Delivery time 1 d

Optional Services

DAkkS Calibration 962-617

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-617	DAkkS Calibration KERN 962-617
962-617R	DAkkS Recalibration KERN 962-617R
920-712	Marking (Last digits of serial number) KERN 920-712

SPARE PARTS

Model	Description
314-040-600	Aluminium Protective Box KERN 314-040-600
315-242	Tweezers KERN 315-242