

PRODUCT BROCHURE

KERN 367-01

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M3, 1 g, Knob, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	367-01
Model-Serie	367-01
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761257207
Predecessor Model	367-41
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g
Tolerance (OIML)	10 mg
Case aluminium	317-010-600
Case plastic	347-030-400
Case wood	337-010-200

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	1 g
---------------	-----

OIML class	M3
OIML conform	yes
Tolerance (OIML)	10 mg
Density (kg/m ³)	7950
Volume (cm ³)	0,119

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	6×6×6 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	no
-----------------------	----

Optional Services

DAkkS Calibration	962-631
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>337-010-200</u>	Wooden Box for Set of Weights KERN 337-010-200
<u>347-030-400</u>	Plastic Weight Box KERN 347-030-400
<u>317-010-600</u>	Aluminium Protective Box KERN 317-010-600
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>962-631</u>	DAkkS Calibration KERN 962-631
<u>962-631R</u>	DAkkS Recalibration KERN 962-631R