

PRODUCT BROCHURE

KERN 347-066

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M1, 50 g, Hook weight, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class M1, 50 g, Hook weight, stainless steelfine turned (OIML)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-066
Model-Serie	347-066
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761256200
Predecessor Model	347-466
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	50 g
Tolerance (OIML)	3 mg
Case plastic	347-070-400

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	50 g
OIML class	M1
OIML conform	yes
Tolerance (OIML)	3 mg
Density (kg/m ³)	7950
Volume (cm ³)	5,9524

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	18×18×36 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-636
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>347-070-400</u>	Plastic Weight Box KERN 347-070-400
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>962-636</u>	DAkkS Calibration KERN 962-636
<u>962-636R</u>	DAkkS Recalibration KERN 962-636R