


PRODUCT BROCHURE

KERN 347-127

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M1, 2 kg, Slot weight (Type B), stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class M1, 2 kg, Slot weight (Type B), stainless steelfine turned (OIML)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-127
Model-Serie	347-127
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366794
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	2 kg
Tolerance (OIML)	100 mg

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	2 kg
OIML class	M1

OIML conform	yes
Tolerance (OIML)	100 mg
Density (kg/m ³)	8000
Volume (cm ³)	238,0952

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	118×118×31 mm
Gross weight	2 kg
Net weight	2 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-642
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-280</u>	Gloves KERN 317-280
<u>347-137-100</u>	Supporting Bar KERN 347-137-100
<u>347-147-100</u>	Supporting Bar KERN 347-147-100
<u>347-117-100</u>	Supporting Bar KERN 347-117-100

SERVICES

Model	Description
<u>962-642R</u>	DAkkS Recalibration KERN 962-642R
<u>962-642</u>	DAkkS Calibration KERN 962-642