

PRODUCT BROCHURE

KERN 347-126-N00

Date: 28.05.2026

ENGLISH



Individual weight (Newton), 20 N, Hook weight, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-126-N00
Model-Serie	347-126-N00
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 147502 1
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Tolerance (OIML)	100 mg
Case plastic	347-130-400

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

OIML class	M1
OIML conform	no
Tolerance (OIML)	100 mg
Density (kg/m ³)	7950
Volume (cm ³)	238,0952

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	60×60×146 mm
Gross weight	2 kg
Net weight	2 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	6 d

Optional Services

DAkkS Calibration	962-642N
-------------------	----------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
317-280	Gloves KERN 317-280
347-130-400	Plastic Weight Box KERN 347-130-400

SERVICES

Model	Description
962-642RN	DAkkS Recalibration KERN 962-642RN
962-642N	DAkkS Calibration KERN 962-642N

ENGLISH

RELATED ARTICLES DATA

Model name	Tolerance (OIML)	Case plastic
347-126-N00	100 mg	347-130-400