

KERN ODC 241

Tablet with integrated camera for the very best monitoring of the sample as well as digital documentation



Article codes, product groups & Status information

Model article-code	ODC 241
Model-Serie	ODC-2
Assortment	Optics
Product Type	Microscopes
Product group	Microscope cameras
Product usage type	Main article
EAN code (Model)	4045761217683
Planned start date	2018-2HJ
Customs tariff number	85258900

Technical data - Primary

Camera type	Tablet camera
Camera resolution	5 MP
Sensor type	Sony CMOS
Sensor size	1/2,8"

Technical data - Construction

Color capability	Colour
Camera mounting	C-Mount
Interfaces	Bluetooth;HDMI;SD Kartenslot;USB 2.0;WLAN
Operating system	Android
Applicable software	Integrated
Software functions	Image capturing;Measurements;Video capturing

Dimensions (W×D×H)	238×182×51 mm
--------------------	---------------

Technical data - Display

Display type	LCD touch
Screen size	9.7"
Display resolution	2048x1536

Technical data - Power supply

Power supply kind	Power adapter
Input Voltage	100-240V AC 50/60Hz 0.5A
Power Plug	Inside plus, outside minus
Protection Class	Class III - Low voltage
Cable attachment	removable

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	275×230×85 mm
Gross weight	1.2 kg
Net weight	0.80 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	206.483 g
Packaging component - by weight - plastics	1.539 g
Packaging component - by weight - styrofoam	48.4785 g

KERN ODC 241

Tablet with integrated camera for the very best monitoring of the sample as well as digital documentation



Packaging component - by weight - foam 0 g

Delivery time 1 d

Technical data - Environmental conditions

Storage Temp Min -5 °C

Storage Temp Max 40 °C

Technical data - Verification

WEEE Yes

Functions

Standard

