

## PRODUCT BROCHURE

---

# KERN 315-311-306

Date: 19.09.2025

ENGLISH



Safety Set 1 kg (F2); 50 g (F2)

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## DESCRIPTION

- These “KERN Safety Sets” each contain one test weight for testing the sensitivity, i.e. the correct adjustment of your balance as well as a small test weight to test the lower end of the weighing range, the so-called minimum sample weight
- As an option, the “KERN Safety Set” has space for another test weight, for testing your balance at a weight which is relevant for you
- Useful accessories which have been selected to suit that particular “KERN Safety Set”, such as, for example, special gloves, tweezers, weight grips, brushes, etc. will assist you in handling your test weights properly
- Stored in the practical protective case next to your balance, you can check and ensure the high precision of your balance at any time

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	315-311-306
Model-Serie	315
Assortment	KERN
Product Type	Test weights
Product group	Safety sets
Product usage type	Main item;Accessories/Services

EAN code (Model) 4045761223424

Customs tariff number 84239090

### Technical data - Measuring system

Nominal weights, included 1 kg (F2); 50 g (F2)

Nominal value 50 g & 1 kg

OIML class F2

OIML conform yes

## Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	255×135×140 mm
Gross weight	3 kg
Net weight	1.2 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	1.8 kg
Delivery time	1 d

## Technical data - Verification

Verification possible	yes
-----------------------	-----

## Optional Services

DAkkS Calibration	962-400-400
-------------------	-------------

## FUNCTIONS

Standard
----------



## SERVICES

Model	Description
952-441	Verification KERN 952-441
921-441	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-441
962-400-400	DAkkS Calibration KERN 962-400-400

## SPARE PARTS

Model	Description
315-242	Tweezers KERN 315-242