

KERN 347-05

Individual weight, OIML Class M1, 20 g, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-05
Model-Serie	347
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761008687
Planned start date	2004-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	20 g
Tolerance (OIML)	2,5 mg
Case aluminium	317-050-600
Case plastic	347-050-400
Case wood	337-050-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	20 g
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	2,5 mg
Density (kg/m ³)	7950

Functions

Standard

Volume (cm³) 2,5

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	13×13×23 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-635
-------------------	---------

KERN 347-05

Individual weight, OIML Class M1, 20 g, Knob, stainless steel fine turned

