



Matrix Needle Printer KERN 911-013*

For printing weighing results on normal paper, which is ideal for long-term archiving.



Matrix Needle Printer KERN YKN-01*

For printing weighing results on normal paper, which is ideal for long-term archiving.



Thermal Printer KERN YKC-01

Reliable thermal printer for standard applications.

Technical data

- Needle print
- Speed: 1,8 line/s
- Max. 40 characters per line
- Width of the paper strip: 57 mm wide
- Operation via mains adapter 230 V/50 Hz (standard EU), as standard
- Dimensions housing: W×D×H 100×180×106 mm
- Not suitable for KERN PNS/PNJ, PES/PEJ ranges
- Interface cable not included with delivery

Price printer without interface cable,
KERN 911-013

Accessories

- Paper rolls (10 pcs.), KERN 911-013-010
- Power Adapter, EU/GB/US/CH/AUS, KERN YKA-42

Technical data

- Needle print
- Speed: 1 line/s
- Max. 16 characters per line
- Width of the paper strip: 44 mm wide
- Operation via mains adapter 230 V/50 Hz (standard EU/GB/US), as standard
- Dimensions housing: W×D×H 158×106×40 mm
- Interface cable not included with delivery

Price printer without interface cable,
KERN YKN-01

Accessories

- Paper rolls (5 pcs.), KERN YKN-A01

Technical data




- Thermal print
- Speed: 300 mm/s
- Max 58 characters per line
- Supported paper widths: 80 mm (roll)
- Operation via mains adapter 230 V/50 Hz standard EU, as standard
- Interfaces: RS-232
- Dimensions housing: W×D×H 125×125×150 mm
- Interface cable not included with delivery

Price printer without interface cable,
KERN YKC-01

Accessories

- Paper rolls (5 pcs.), KERN YKC-A10

! * ONLY WHILE STOCKS LAST



Important note for printer orders:

Please state your balance model when ordering, so that the printer can be preconfigured in the factory to suit your balance.

Individual header data:

For printers **KERN 911-013** and **YKN-01** the free software **KERN SHM-01** (see www.kern-sohn.com) can be used to print 4 header lines on the printout (see marking on balance).
Compatible models: KERN 440, DE, FCB, PCD.