


KERN 347-127-N00

Individual weight (Newton), 20 N, Slot weight (Type B), stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-127-N00
Model-Serie	347-7-N00
Assortment	Scales
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366800
Planned start date	2024-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Tolerance (OIML)	100 mg
------------------	--------

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

OIML class	M1
Density (kg/m ³)	8.000
Volume (cm ³)	238,0952

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	118×118×31 mm
Gross weight	2,0 kg
Net weight	2,0 kg
Shipping method	Parcel service

Functions

Packaging component - by weight - plastics 0 g

Delivery time 6 d

Optional Services

DAkkS Calibration 962-642N

KERN 347-127-N00

Individual weight (Newton), 20 N, Slot weight (Type B), stainless steel fine turned



Standard



Option

