

PRODUCT BROCHURE

KERN 344-066

Date: 10.03.2026

ENGLISH



Weight set, M1, 1 g - 1 kg, Knob, stainless steel fine turned, in aluminium case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	344-066
Model-Serie	344
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	404576 122484 1
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 1 kg
Case aluminium	314-060-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Number of weights in set	13
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 g - 1 kg
OIML class	M1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	230×140×110 mm
Gross weight	3.6 kg
Net weight	3.2 kg

Shipping method	Parcel service
Packaging component - by weight - plastics	350 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-619
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-619</u>	DAkkS Calibration KERN 962-619
<u>962-619R</u>	DAkkS Recalibration KERN 962-619R
<u>920-714</u>	Marking (Last digits of serial number) KERN 920-714

SPARE PARTS

Model	Description
<u>314-060-600</u>	Aluminium Protective Box KERN 314-060-600
<u>315-242</u>	Tweezers KERN 315-242