

PRODUCT BROCHURE

KERN 348-22

Date: 11.03.2026

ENGLISH



Set of milligram weights, M1, 1 mg - 500 mg, Platelet, stainless steel polished, in plastic case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	348-22
Model-Serie	348
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1495074
Predecessor Model	348-22-2024a
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	1 mg - 500 mg
Case plastic	328-222-420
Weight set case	in plastic case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	12

Technical data - Measuring system

OIML class	M1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	136×53×22 mm
Gross weight	0.10 kg
Net weight	0.10 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	0 g

Delivery time 1 d

Optional Services

DAkkS Calibration 962-650

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>338-090-200</u>	Wooden Box for Set of Weights KERN 338-090-200
<u>347-009-400</u>	Plastic Weight Box KERN 347-009-400
<u>317-009-600</u>	Aluminium Protective Box KERN 317-009-600
<u>348-01</u>	Test Weight KERN 348-01
<u>348-02</u>	Test Weight KERN 348-02
<u>348-03</u>	Test Weight KERN 348-03
<u>348-04</u>	Test Weight KERN 348-04
<u>348-05</u>	Test Weight KERN 348-05
<u>348-06</u>	Test Weight KERN 348-06
<u>348-07</u>	Test Weight KERN 348-07
<u>348-08</u>	Test Weight KERN 348-08
<u>348-09</u>	Test Weight KERN 348-09
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270

SERVICES

Model	Description
<u>962-650</u>	DAkkS Calibration KERN 962-650
<u>962-650R</u>	DAkkS Recalibration KERN 962-650R
<u>920-701</u>	Marking (Last digits of serial number) KERN 920-701

SPARE PARTS

Model	Description
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270