

PRODUCT BROCHURE

KERN PCD-A04

Date: 13.12.2025

ENGLISH



Internal rechargeable battery pack, operating time up to 24 h, charging time approx. 12 h, AUTO-OFF function to preserve the rechargeable battery pack

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	PCD-A04
Model-Serie	PCD
Assortment	KERN
Product Type	Laboratory balances
Product group	Accessories (printer, etc.)
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761134188
Customs tariff number	85075000
CE Label	yes

Technical data - Power supply

Mains adapter	Battery & Power supply
---------------	------------------------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	110×85×45 mm
Gross weight	0.25 kg
Net weight	0.20 kg
Shipping method	Parcel service
UN number	3496
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	20 g
Packaging component - by weight - plastics	20 g

Packaging component - by weight - styrofoam 0 g

Packaging component - by weight - glass 0 g

Packaging component - by weight - wood 0 g

Packaging component - by weight - foam 0 g

Packaging component - by weight - metal 0 g

Delivery time 1 d

FUNCTIONS

Standard

