

KERN 326-09

Individual weight, OIML Class F1, 500 g, Eco Shape, stainless steel polished



Article codes, product groups & Status information

Model article-code	326-09
Model-Serie	326
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761223660
Predecessor Model	326-09-2018a; 327-69; 329-09
Successor Model	326-09-2025e
Planned start date	2019-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	500 g
Tolerance (OIML)	2,5 mg
Case aluminium	317-090-600
Case plastic	347-090-400
Case wood	317-090-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	500 g
OIML class	F1
OIML conform	Yes

Functions

Tolerance (OIML)	2,5 mg
------------------	--------

Density (kg/m ³)	7950
------------------------------	------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	38×38×64 mm
------------------------------	-------------

Gross weight	0.50 kg
--------------	---------

Net weight	0.50 kg
------------	---------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - plastics	0 g
--	-----

Delivery time	1 d
---------------	-----

Optional Services

DAkKS Calibration	962-439
-------------------	---------

KERN 326-09

Individual weight, OIML Class F1, 500 g, Eco Shape, stainless steel polished



Standard

