

## PRODUCT BROCHURE

---

# KERN 317-16

Date: 09.03.2026

ENGLISH



Individual weight, OIML Class E2, 50 kg, Knob, stainless steel polished

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

- Individual weight, OIML Class E2, 50 kg, Knob, stainless steel polished (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section
- For the verification you need a box. Please add a box for single weights.

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 317-16                         |
| Model-Serie           | 317                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML E2 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 404576 1003828                 |
| Customs tariff number | 84239090                       |
| Country of Origin     | DE                             |

### Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 50 kg       |
| Tolerance (OIML) | 80 mg       |
| Case aluminium   | 317-160-600 |
| Case wood        | 317-160-100 |

### Technical data - Construction

|          |                          |
|----------|--------------------------|
| Material | stainless steel polished |
|----------|--------------------------|

### Technical data - Measuring system

|                              |       |
|------------------------------|-------|
| Nominal value                | 50 kg |
| OIML class                   | E2    |
| OIML conform                 | yes   |
| Tolerance (OIML)             | 80 mg |
| Density (kg/m <sup>3</sup> ) | 7950  |
| Volume (cm <sup>3</sup> )    | 6250  |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 178×178×298 mm |
| Gross weight                               | 50 kg          |
| Net weight                                 | 50 kg          |
| Shipping method                            | Freight        |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 30 d           |

### Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-346 |
|-------------------|---------|

## FUNCTIONS

Standard



Option



## ACCESSORIES

| Model                              | Description                                    |
|------------------------------------|------------------------------------------------|
| <a href="#"><u>317-290</u></a>     | Gloves KERN 317-290                            |
| <a href="#"><u>317-160-100</u></a> | Wooden Box for Set of Weights KERN 317-160-100 |
| <a href="#"><u>317-160-600</u></a> | Aluminium Protective Box KERN 317-160-600      |

## SERVICES

| Model                           | Description                                          |
|---------------------------------|------------------------------------------------------|
| <a href="#"><u>962-346</u></a>  | DAkkS Calibration KERN 962-346                       |
| <a href="#"><u>962-346R</u></a> | DAkkS Recalibration KERN 962-346R                    |
| <a href="#"><u>921-346</u></a>  | Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-346 |