

KERN 337-141

Individual weight, OIML Class F2, 10 kg, Check weight, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	337-141
Model-Serie	337
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761006577
Planned start date	2004-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	160 mg
Case aluminium	Not available
Case plastic	Not available
Case wood	337-141-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	10 kg
OIML class	F2
OIML conform	Yes
Tolerance (OIML)	160 mg
Density (kg/m ³)	7950

Functions

Standard

Volume (cm³) 1250

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	136×136×132 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-444
-------------------	---------

KERN 337-141

Individual weight, OIML Class F2, 10 kg, Check weight, stainless steel fine turned

