

KERN 315-312-308

Safety Set 200 g (F2); 2000 g (F2)



Article codes, product groups & Status information

Model article-code	315-312-308
Model-Serie	315
Assortment	KERN
Product Type	Test weights
Product group	Safety sets
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761223448
Planned start date	2019-1HJ
Customs tariff number	84239090

Technical data - Measuring system

Nominal weights, included	200 g (F2); 2000 g (F2)
Nominal value	200 g & 2 kg
OIML class	F2
OIML conform	Yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	255×135×140 mm
Gross weight	3 kg
Net weight	1.2 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	1.8 kg
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-400-400
-------------------	-------------

Functions

KERN 315-312-308

Safety Set 200 g (F2); 2000 g (F2)



Standard

