

KERN ODC 241

Tablet with integrated camera for the very best monitoring of the sample as well as digital documentation



Article codes, product groups & Status information

Model article-code	ODC 241
Model-Serie	ODC-2
Assortment	Optics
Product Type	Microscopes
Product group	Microskope cameras
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761217683
Planned start date	2018-2HJ
Customs tariff number	85258900
CE Label	Yes

Technical data - Primary

Camera type	Tablet camera
Camera resolution	5 MP
Sensor type	Sony CMOS
Sensor size	1/2,8"

Technical data - Construction

Color capability	Colour
Camera mounting	C-Mount
Interfaces	Drahtlose Datenübertragung;HDMI;HDMI Micro;SD Kartenslot;USB 2.0;WLAN
Operating system	Win 10;Win 11
Applicable software	Integrierte Kamerasoftware

Software functions

Image capturing;Measurements;Video capturing

Dimensions (W×D×H)

238×182×51 mm

Technical data - Display

Display type	LCD touch
Screen size	9,7"
Display resolution	2.048x1.536

Technical data - Power supply

Power supply kind	Power adapter
Input Voltage	100-240V AC 50/60Hz 5A
Power Plug	Inside plus, outside minus
Protection Class	Class III - Low voltage
Cable attachment	Removable

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	275×230×85 mm
Gross weight	1,2 kg
Net weight	0,80 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	206,483 g
Packaging component - by weight - plastics	1,539 g

KERN ODC 241

Tablet with integrated camera for the very best monitoring of the sample as well as digital documentation



Packaging component - by weight - styrofoam	48,4785 g
---	-----------

Packaging component - by weight - foam	0 g
--	-----

Delivery time	1 d
---------------	-----

Technical data - Environmental conditions

Storage Temp Min	-5 °C
------------------	-------

Storage Temp Max	40 °C
------------------	-------

Technical data - Verification

WEEE	Yes
------	-----

Functions

Standard

