

PRODUCT BROCHURE

KERN 303-05

Date: 09.03.2026

ENGLISH



Weight set, E1, 1 mg - 500 g, Knob, stainless steel polished, in wooden case


CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information



ENGLISH

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

Model article-code	303-05
Model-Serie	303
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main article
EAN code (Model)	4045761002494
Customs tariff number	84239090
Country of Origin	DE

Nominal value	1 mg - 500 g
---------------	--------------

Case wood	313-050-100
Weight set case	in wooden case

Material	stainless steel polished
Is weight set	yes
Number of weights in set	24
Type of case	in wooden case

Nominal value	1 mg - 500 g
OIML class	E1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	230×180×110 mm
Gross weight	3.2 kg
Net weight	3.2 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-204
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>963-204</u>	DAkkS Calibration KERN 963-204
<u>962-204R</u>	DAkkS Recalibration KERN 962-204R
<u>920-713</u>	Marking (Last digits of serial number) KERN 920-713

SPARE PARTS

Model	Description
<u>313-050-100</u>	Box of Weights KERN 313-050-100
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270