

PRODUCT BROCHURE

KERN

3 15-3 12-407

Date: 19.03.2025

ENGLISH



Safety Set 100 g (M1); 2000 g (F2)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

- These “KERN Safety Sets” each contain one test weight for testing the sensitivity, i.e. the correct adjustment of your balance as well as a small test weight to test the lower end of the weighing range, the so-called minimum sample weight
- As an option, the “KERN Safety Set” has space for another test weight, for testing your balance at a weight which is relevant for you
- Useful accessories which have been selected to suit that particular “KERN Safety Set”, such as, for example, special gloves, tweezers, weight grips, brushes, etc. will assist you in handling your test weights properly
- Stored in the practical protective case next to your balance, you can check and ensure the high precision of your balance at any time

SPECIFICATIONS

Article codes, product groups & Status information

| | |
|--------------------|--------------------------------|
| Model article-code | 315-312-407 |
| Model-Serie | 315 |
| Assortment | KERN |
| Product Type | Test weights |
| Product group | Safety sets |
| Product usage type | Main item;Accessories/Services |

EAN code (Model) 4045761245747

Customs tariff number 84239090

Technical data - Measuring system

Nominal weights, included 100 g (M1); 2000 g (F2)

Nominal value 100 g & 2 kg

OIML class F2;M1

Technical data - Packing & Shipping

| | |
|--|----------------|
| Dimensions packaging (W×D×H) | 255×135×140 mm |
| Gross weight | 0.05 kg |
| Net weight | 0.05 kg |
| Shipping method | Parcel service |
| Packaging component - by weight - plastics | 0 g |
| Delivery time | 1 d |

Technical data - Verification

| | |
|-----------------------|-----|
| Verification possible | yes |
|-----------------------|-----|

Optional Services

| | |
|-------------------|-------------|
| DAkkS Calibration | 962-400-600 |
|-------------------|-------------|

FUNCTIONS

| |
|----------|
| Standard |
|----------|



SERVICES

| Model | Description |
|-------------|------------------------------------|
| 962-400-600 | DAkkS Calibration KERN 962-400-600 |

SPARE PARTS

| Model | Description |
|---------|-----------------------|
| 315-242 | Tweezers KERN 315-242 |