

KERN YKA-12

Date: 28.01.2025

ENGLISH



Power adapter (EU/UK/US)

CONTACT

 Address

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0
Fax : +49 - [0] 7433 - 9933 - 149

 Internet

Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Power adapter YKA-12 (9 V, 800 mA | 100 V - 240 V, 50 / 60 Hz | EURO, UK, US, CH)

SPECIFICATION

Article codes, product groups & Status information

Model article-code	YKA-12
Model-Serie	YKA
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc)
Product usage type	Spare part
EAN code (Model)	4045761178373
Planned start date	2016-1HJ
Customs tariff number	85044095

Packaging component -
by weight - glass 0 g

Packaging component -
by weight - wood 0 g

Delivery time 1 d

ENGLISH

Technical data - Power supply

Mains adapter/adaptor included	CH;EURO;UK;US
Input Voltage	100 V - 240 V, 50 / 60 Hz

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	90×80×50 mm
Gross weight	0.40 kg
Net weight	0.35 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	15 g
Packaging component - by weight - plastics	15 g

RELATED ARTICLES DATA

Model name	Max. Force
YKA-12	
824-902	
RYKB-5001	
RYDB-9003	
RTYPS-8001	
RTYPS-9001	
RTYPS-8002	
RTMCN-9001	
YMM-05	
RYKB-8001	
RYKS-8001	
RYKU-5001	
RYKU-5002	
RYKN-5001	