

PRODUCT BROCHURE

OPTICS OBB-A1227

Date: 20.03.2025

ENGLISH



Monocular lens barrel, 30° angled, 360° revolving

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Microscope tube monocular OBB-A1227

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1227
Model-Serie	OBB-AT
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761167599
Customs tariff number	90119000

Dimensions packaging (W×D×H)	190×130×100 mm
------------------------------	----------------

Gross weight	1.8 kg
--------------	--------

Net weight	1.2 kg
------------	--------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - plastics	600 g
--	-------

Delivery time	1 d
---------------	-----

FUNCTIONS

Standard



ENGLISH

Technical data - Primary

Tube	Monocular
------	-----------

Technical data - Construction

Tube angle	30°
Tube rotation	yes

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping