

PRODUCT BROCHURE

OPTICS OBB-A2521

Date: 16.03.2025

ENGLISH



Object holder for 96-hole microtiter plate for inverted microscope

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Object holder for 96-hole microtitre plate OBB-A2521

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A2521
Model-Serie	OBB-AZ
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761173453
Customs tariff number	90119000

Shipping method Parcel service

Packaging component - by weight - plastics 200 g

Delivery time 1 d

FUNCTIONS

Standard



ENGLISH

Technical data - Construction

Dimensions	150×100×5 mm
------------	--------------

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	145×120×10 mm
Gross weight	0.30 kg
Net weight	0.10 kg