

PRODUCT BROCHURE

**OPTICS
OBB-A1564**

Date: 16.03.2025

ENGLISH

Achromatic objective, E-Plan 40 x /0,65 (sprung) W.D. ((0,85 mm), anti-fungus

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Achromatic objective, E-Plan 40 x /0,65 (sprung) W.D. (0,85 mm), anti-fungus OBB-A1564

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1564
Model-Serie	OBB-AO
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761229013
Customs tariff number	90119000

FUNCTIONS

Standard



Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	65×40×40 mm
Gross weight	0.60 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	300 g
Delivery time	1 d