

PRODUCT BROCHURE

KERN 346-84-2024e

Date: 19.03.2025

ENGLISH



Image not
available

Individual weight, OIML Class M1, 1000 kg, Heavy Duty Weight, cast iron lacquered

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- Individual weight, OIML Class M1, 1000 kg, Heavy Duty Weight, cast iron lacquered
- For the verification you need a box. Please add a box for single weights.

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	346-84-2024e
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 test weights
Product usage type	Hauptartikel;Zubehör/DL
EAN code (Model)	404576 1479586
Predecessor Model	346-84
Planned start date	2024-2HJ

Technical data - Primary

Nominal value	1000 kg
Tolerance (OIML)	50 g
Case aluminium	nicht verfügbar
Case plastic	not available
Case wood	nicht verfügbar

Technical data - Construction

Dimensions	750×500×500 mm
Material	cast iron lacquered

Technical data - Measuring system

Nominal value	1000 kg
OIML class	M1
OIML conform	yes
Tolerance (OIML)	50 g
Density (kg/m ³)	7100
Volume (cm ³)	138889

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	750×500×500 mm
Gross weight	1000 kg
Net weight	1000 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	2 d

Optional Services

DAkKS Calibration	962-694
-------------------	---------

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium
346-06	5 kg	250 mg	346-060-600
346-07	10 kg	500 mg	346-070-600
346-08	20 kg	1 g	346-080-600
346-09	50 kg	2,5 g	346-090-600
346-86	5 kg	250 mg	346-060-600
346-87	10 kg	500 mg	346-070-600
346-77	10 kg	500 mg	346-070-600
346-76	5 kg	250 mg	346-060-600
346-79-2023a	50 kg	2,5 g	-
346-88	20 kg	1 g	346-080-600
346-78	20 kg	1 g	346-080-600
346-89	50 kg	2,5 g	-
346-81	100 kg	5 g	-
346-82	200 kg	10 g	-
346-83	500 kg	25 g	-
346-84	1000 kg	50 g	-
346-022-005	20 kg	1 g	-
346-042-008	40 kg	1,5 g	-
346-052-010	50 kg	2,5 g	-
346-055-004	50 kg	2,5 g	-
346-055-009	50 kg	2,5 g	-
346-065-009	60 kg	3 g	-
346-79	50 kg	2,5 g	346-090-600
346-86-2025e	5 kg	250 mg	346-060-600
346-87-2025e	10 kg	500 mg	346-070-600
346-88-2025e	20 kg	1 g	346-080-600
346-89-2025e	50 kg	2,5 g	346-090-600